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## HE staining experiment report

### 1. Laboratory equipment and reagents

#### 1.1 Laboratory equipment

Equipment	Manufacturers	Model
Dehydrator	Hubei Benor Medical Technology Co., Ltd	JT-12S
Automatic biological tissue embedding machine	Hubei Benor Medical Technology Co., Ltd	JB-L5
Paraffin embedding machine (Cold Plate )	Hubei Benor Medical Technology Co., Ltd	JB-P5
Tissue spreader	KEDEE	51A011
Rotary microtome	Shanghai Leica Instrument Co., Ltd	HistoCoreBIO CUT
Oven	Shanghai Yuejin Medical Equipment Co., LTD	HGPF-50
Adhesive slide and Cover glass	Jiangsu Huida medical Equipment Co., LTD	710510
Dyeing machine	DAKEWE	DP260s
Upright optical microscope	NIKON	ECLIPSE E100

#### 1.2 Laboratory reagents

Reagents	Manufacturers	Catlog
Anhydrous ethanol	Hangzhou Hongda Chemical Instrument	SJ003614

---

	Co., LTD	
Xylene	Sinopharm Chemical Reagent co., Ltd.	10023418
Hematoxylin dye	Hangzhou Haoke Biotechnology Co., Ltd	HK2053
Acid Alcohol		HK2054
Differentiation Solution	Hangzhou Haoke Biotechnology Co., Ltd	
Eosin	Dakewe	HK2051
Ammonia Solution	Hangzhou Haoke Biotechnology Co., Ltd	HK2055
Neutral gum	Sinopharm Chemical Reagent co., Ltd.	10004160

---

## 2. Dyeing steps

**2.1 Dewaxing and hydration:** Put the sections into xylene I 8 min - xylene II 8 min - xylene III 8 min - anhydrous ethanol I 5 min - anhydrous ethanol II 5min - 85% Ethyl alcohol for 5 min - 75% Ethyl alcohol for 5 min, wash with tap water.

**2.2 Hematoxylin staining:** Put sections into Hematoxylin solution for 6 min. Then treat the section with acid alcohol differentiation solution for 2 s, rinse with tap water. Treat the section with Hematoxylin Bluing solution for 15-30 s, rinse with tap water.

**2.3 Eosin staining:** Place the sections in 95% ethanol for 1 min, eosin dye for 10-30s.

**2.4 Dehydration and sealing:** Put the sections into absolute ethanol I for 30 s - absolute ethanol II for 2.5 min - absolute ethanol III for 2.5 min - xylene I for 2.5 min - xylene II for 2.5 min ,sealing with neutral gum.

**2.5** Microscope inspection, image acquisition and analysis.

## 3. Interpretation of results:

The nucleus is blue and the cytoplasm is red.

#### **4. Precautions:**

4.1 Pay attention to whether the section is thoroughly dewaxed;

4.2 Pay attention to the differentiation degree of nucleus;

4.3 Pay attention to the potency of hematoxylin and eosin, and change the dyeing solution in time.